
## 1. GERAL

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## 1. Overview

- Product Category : Glass-ceramic Ingot
- Product Name : Rosetta SP
- Intended Use of the Device : This product is used to make artificial tooth in order for damaged tooth to get its function and beauty restored.
- Packaging Unit : Refer to HASS standard package.

## 2. Instruction

Ingot	Wax + Sprue	Invest Ring	**
R10 1ea	up to 0,7 g	100 g	** Refer to the pressing schedule
R20 1ea	up to 1,7 g	200 g	*** Do not use the alumina beads for the sand blast, just use the glass beads

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## 3. Cautions

### (1) Cautions before Use

- ① Check the defects of the product for any damage or crack before use.
- ② Take care of burn when inserting Ingot into investment ring.
- ③ Make sure plunger is well applied and dried up with parting agent before inserting.
- ④ Inserting Ingot / plunger, and investment into to press furnace should be finished in a short time.
- ⑤ Cool down up to the room temperature after completing pressing procedure.
- ⑥ Be careful not to inhale dust during divesting process.
- ⑦ Handle with care not to damage the pressed product when using wheel and bur.

### (2) Storage and Maintenance before Use

- ① Store the product at room temperature in a dry place.
- ② Pack and store the product properly to ensure that it is not damaged.
- ③ Store the product at temperatures ranging from 0°C~40°C, in combination with relative humidity of 10% r.H ~ 90% r.H, under atmospheric pressures ranging from 500 hPa ~ 1060 hPa.

## 4. Mechanical and Physical Properties

- ① Material : Glass-ceramics
- ② Flexural Strength : over 300 MPa
- ③ Chemical Solubility : below 100 µg/cm<sup>2</sup>
- ④ Coefficient of Thermal Expansion : 10.1 (±0.5) x 10<sup>-6</sup> K<sup>-1</sup>

\* This is a single-use product. \* Do not reuse.

## 5. Pictograph

do not reuse	Caution	Catalogue Number	Consult instructions for Use
Batch Code	Date of Manufacture	Do not use if package is damaged	
Manufacturer	Authorized Representative in the European Community	Non Sterile	

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## ⓘ Caution

There may be a slight difference between the displayed temperature and the actual temperature of each furnace. Before you inject the above recommended schedule is suitable for the furnace being used. Otherwise, try to find the optimized pressing temperature through the following process.

1) If there are some traces of tiny bubbles on the surface of the pressed restoration, reduce the maximum temperature by 5-10°C and retry the pressing procedure.

2) If the marginal area of restoration is not formed completely, increase the maximum temperature by 5-10°C and retry the pressing procedure.

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